



KEY FEATURES

- Castellated internal foam provides flexibility when used on a profiling bed
- Hygienic wipe-clean surface
- Water-resistant welded seams
- Anti tamper properties
- Reversible - both sides
- CMHR soft foam
- Compliant with BS 7177:2008 + A1:2011 Crib Source 5



MATERIAL/FABRICATION

Designed for compatibility with a profiling bed. The mattresses internal foam castellations allow flexibility. PU Coated modacrylic cotton for comfort and CMHR foam for safety. Water-resistant, HF welded seams to prevent liquid penetration.

TEST STANDARDS

Fire retardant to BS EN 597-1-2:1995, BS 6807:2006 and BS 7177 2008 + A1 2011 Crib Source 5 for medium hazard use.

OPTIONS

- Solo Profiling beds available (suitable for use with Flex Profile Mattress)

COLOUR/FINISH

Cover is Navy Blue as standard

INSTALLATION/FITTING

No fitting required

MAINTENANCE

Wipe with mild detergent. Wash or dry-clean using detergent up to 40°C. For disinfection use sodium hypochlorite up to concentrations of 10,000 parts per million of active chlorine. If the fabric has been heavily soiled a chlorine solution of up to 10,000 ppm can be used. Frequent or prolonged use of 1% chlorine solution (10,000 ppm) may reduce the useful life of the fabric. Avoid using biological or phenolic-based cleaners. Wash and dry completely before storage.

WARRANTY

2 year guarantee

